DE

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HINWEISE

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OPERATING MANUAL

NOTE

- Please read the user guide before using the unit.
- Do not place the device in water – it is not waterproof.
- Do not clean the tip with hot water; it could get damaged.
- After each reading clean the tip carefully with a dry cloth and disinfect with ethyl alcohol to prevent incorrect readings being taken.
- Never expose the BS 31 to direct sunlight or extreme heat, as this might damage the display.
- Never attempt to bend the thermometer using force.
- Only open the device when the batteries need changing.
- On rectal temperature readings the device must be inserted carefully and no further than 1.5 cm. There is a risk of injury!
- Do not leave in the vicinity of children or animals.

Normal body temperature lies between 35.8 °C and 38.0 °C. Ideally you should take your own temperature and that of your family when you are fully fit, and write down the reading, so that in the event of illness you will be able to determine straightaway whether the patient has a temperature or not.

OPERATION

1. Press the on/off switch. You will hear a beep.
2. The last temperature reading appears for about 2 seconds on the LCD display.
3. The display then flashes as it prepares the device for instruction manual.
4. Once you hear a beep and the Lo °C display appears continuously, the device is ready for use.

5. Now the measuring tip can be used to take oral or rectal readings, or readings from under the armpits.

6. It takes between 30 and 60 seconds to take a reading. The end of the reading is signalled by a continuous 10-second beep.

Beeps:

– measured temperature < 37.8 °C:
10 long beeps

– measured temperature ≥ 37.8 °C:
30 short beeps

7. To take another reading, switch the BS 31 off for a moment and then on again.

The device switches off automatically 10 minutes after the last reading. However, to save the battery, always use the on/off switch to switch the device off after use.

CHANGING THE BATTERY

The battery must be changed as soon as the battery symbol appears on the display, or the device no longer switches on.

To change the battery, simply remove the end cap on the lid of the battery compartment.

Now you can take out the discharged battery, insert a new one (type LR41/AG3/UCC392/SR41) ensuring correct polarity, and put the cover back on.

NOTES ON DISPOSAL

Old units, marked with the symbol as illustrated, may not be disposed of in the household rubbish.

CARE AND WARRANTY

Before cleaning the unit, disconnect it if necessary from other components; do not use aggressive cleaning agents.

The unit has been carefully checked for defects. If nevertheless you do have cause for complaint, please send us the unit with your proof of purchase. We offer a 2 year warranty from date of purchase.

We are not liable for damage arising from incorrect handling, improper use or wear and tear.

We reserve the right to make technical modifications.

Dear customer, In the case of a defect, please return the device together with the receipt and original packing material to the point-of-sale.

TECHNICAL DATA

Temperature range: 32 °C-42 °C

Display: degrees celsius

Memory: last reading automatically stored

Type of battery: LR41, AG3, UCC392 or SR41

Battery life: 2 years or 2500 measurements

Automatic cut-out: the device switches off automatically 10 minutes after the reading is taken.

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